

Isobar Vignette - iTRAQ and TMT Analysis

Florian P. Breitwieser, Jacques Colinge

October 7, 2011

Contents

1	Introduction	1
2	Loading data	1
2.1	ibspiked test samples	2
2.2	MSnbase integration	3
3	Data Analysis	3
3.1	Reporter mass precision	3
3.2	Normalization and isotope impurity correction	3
3.3	Fitting a noise model	6
3.4	Protein and peptide ratio calculation	6
3.5	Protein ratio distribution and selection	9
3.6	Detection of proteins with no specific peptides	11
4	Report generation	11
A	File formats	13
A.1	ID CSV file format	13
A.2	IBSpectra CSV file format	13
B	properties.conf for report generation	13
C	Dependencies	16
C.1	L ^A T _E X and PGF/TikZ	16
C.2	Perl	16
D	Session Information	17

1 Introduction

The **isobar** package is designed as an extensible and interactive environment for data analysis and exploration of data produced by Mass Spectrometry analysis of proteins and peptides labelled with isobaric tags, such as iTRAQ and TMT. **isobar** implements the theory presented in Breitwieser et al., Journal of Proteome Research 2011.

`isobar` allows analyzing iTRAQ 4plex and 8plex, and TMT 2plex and 6plex experiments representing them as `IBSpectra` objects. The respective classes are `iTRAQ4plexSpectra`, `iTRAQ8plexSpectra`, `TMT2plexSpectra`, and `TMT6plexSpectra`.

The first thing you need to do is load the package.

```
R> library(isobar) ## load the isobar package
```

2 Loading data

`isobar` can read identifications and quantifications from tab-separated and MGF files. Perl scripts are supplied to generate a tab-separated version from the vendor formats of Mascot and Phenyx, see appendix C. The format is simple and described in appendix A. Experimental support for the mzIdentML format within R is also available - please contact the maintainer in case of problems.

ID.CSV tab-separated file containing peptide-spectra matches and spectrum meta-information such as retention time, m/z and charge. Generated by parser scripts.

MGF contains peak lists from which quantitative information on reporter tags are extracted. Must be centroided.

IBSPECTRA.CSV tab-separated file containing the same columns as ID.CSV plus *quantitative information* extracted from MGF file - that means the reporter tag masses and intensities as additional columns.

`readIBSpectra` is the primary function to generate a `IBSpectra` object. The first argument is one of `iTRAQ4plexSpectra`, `iTRAQ8plexSpectra`, `TMT2plexSpectra` and `TMT6plexSpectra` and denotes the tag type and therefore class.

```
R> ## generating IBSpectra object from ID.CSV and MGF
R> ib <- readIBSpectra("iTRAQ4plexSpectra",list.files(pattern=".id.csv"),
  list.files(pattern=".mgf"))
R> ## write in tabular IBSPECTRA.CSV format to file
R> write.table(as.data.frame(ib),sep="\t",row.names=F,
  file="myexperiment.ibspectra.csv")
R> ## generate from saved IBSPECTRA.CSV - MGF does not have to be supplied
R> ib.2 <- readIBSpectra("iTRAQ4plexSpectra","myexperiment.ibspectra.csv")
```

In case the MGF file is very big, it can be advantageous to generate a smaller version containing only meta- and quantitative information before import in R. On Linux, the tool `grep` is readily available.

```
egrep '^[A-Z]|^1[12][0-9]\.' BIG.mgf > SMALL.mgf
```

2.1 ibspiked test samples

The examples presented are based on the dataset `ibspiked_set1` which has been designed to test `isobar`'s functionality and searched against the Swissprot human database with Mascot and Phenyx. `ibspiked_set1` is an iTRAQ 4-plex data set comprised of a complex background (albumin- and IgG-depleted human plasma) and spiked proteins. MS analysis was performed in ThermoFisher Scientific LTQ Orbitrap HCD instrument with 2D shotgun peptide separation (see original paper for more details). The samples used for each iTRAQ channel are as follows:

- Depleted human plasma background (>150 protein detected);
- Spiked-in proteins with the following ratios
 - CERU_HUMAN (P00450) at concentrations 1 : 1 : 1 : 1;
 - CERU_RAT (P13635) at concentrations 1 : 2 : 5 : 10;
 - CERU_MOUSE (Q61147) at concentrations 10 : 5 : 2 : 1.

A second data set with ratios 1:10:50:100 is available as `ibspiked_set2` from <http://bininformatics.cemm.oeaw.ac.at/isobar>.

The Ceruplasmins have been selected as the share peptides. Hereafter, we load the data package and the ceru protein IDs are identified via the `protein.g` function, which provides a mean to retrieve data from `IBSpectra` objects.

```
R> data(ibspiked_set1)
R> ceru.human <- protein.g(proteinGroup(ibspiked_set1),"CERU_HUMAN")
R> ceru.rat <- protein.g(proteinGroup(ibspiked_set1),"CERU_RAT")
R> ceru.mouse <- protein.g(proteinGroup(ibspiked_set1),"CERU_MOUSE")
R> ceru.proteins <- c(ceru.human,ceru.rat,ceru.mouse)
```

2.2 MSnbase integration

`MSnbase` by Laurent Gatto provides data manipulation and processing methods for MS-based proteomics data. It provides import, representation and analysis of raw MS data stored in `mzXML`, `mzML` and `mzData` using the `mzR` package and centroided and un-centroided MGF peak lists. It allows to use and preprocess raw data whereas `isobar` requires centroided peak lists. In the future, the `isobar` class `IBSpectra` might be based on or replaced by `MSnbase`'s class `MSnSet`. For now, methods for coercion are implemented:

```
R> as(ibspectra,"MSnSet")
R> as(msnset,"IBSpectra")
```

3 Data Analysis

3.1 Reporter mass precision

The distribution of observed masses from the reporter tags can be used to visualize the precision of the MS setup on the fragment level and used to set the correct window for isolation.

The expected masses of the reporter tags are in the slot `reporterMasses` of the implementations of the `IBSpectra` class. The experimental masses are in the matrix `mass` of `AssayData`; they can also be accessed by the method `reporterMasses(x)`.

```
R> sprintf("%.4f",ibspiked_set1@reporterMasses) ## expected masses

[1] "114.1112" "115.1083" "116.1116" "117.1150"

R> mass <- reporterMasses(ibspiked_set1) ## observed masses
R> apply(mass,2,function(x) sprintf("%.4f",quantile(x,na.rm=TRUE,probs=c(0.025,0.975))))

      114      115      116      117
[1,] "114.1110" "115.1081" "116.1115" "117.1148"
[2,] "114.1116" "115.1087" "116.1120" "117.1153"
```

`reporterMassPrecision` provides a plot of the distribution.

3.2 Normalization and isotope impurity correction

Isotope impurity correction factors are supplied by labelling reagent manufacturers. Default values that can be modified by the user are available in `isobar` and corrections are obtained by simple linear algebra.

Due to differences between samples it is advisable to normalize data before further processing. By default, `normalize` corrects by a factor such that the median intensities all the reporter channels are equal.

See figure 2.

```
R> ib.old <- ibspiked_set1
R> ibspiked_set1 <- correctIsotopeImpurities(ibspiked_set1)
R> ibspiked_set1 <- normalize(ibspiked_set1)
```

3.3 Fitting a noise model

A noise model is a approximation of the expected technical variation based on signal intensity. It is stable for a certain experimental setup and thus can be learned once. Noise is observed directly when comparing identical samples in multiple channels (1:1 iTRAQ/TMT sample) and we can use `ibspiked_set1` background proteins as a 1:1 sample. Therefore we exclude the ceruplasmins before fitting a noise model using `NoiseModel`. See figure 3.

```
R> ib.background <- subsetIBSpectra(ibspiked_set1,protein=ceru.proteins,"exclude")
R> noise.model <- NoiseModel(ib.background)

[1] 0.03423 12.14500 1.43708
```

```
R> print(reporterMassPrecision(ibspiked_set1))
```




Figure 1: Reporter mass precision plot.

```

R> png("fig_maplot_normalize.png",width=6.6,height=3,units="in",res=600,pointsize=8)
R> par(mfrow=c(1,2))
R> maplot(ib.old,channel1="114",channel2="117",ylim=c(0.5,2),
          main="before normalization")
R> abline(h=1,col="red",lwd=2)
R> maplot(ibspiked_set1,channel1="114",channel2="117",ylim=c(0.5,2),
          main="after normalization")
R> abline(h=1,col="red",lwd=2)
R> dev.off()

```

pdf

2




Figure 2: Ratio versus intensity plots ('MA plots') before and after applying normalization.

Though only recommended when sufficient data are available, a method exist for the estimation of a noise model without a 1:1 dataset. It takes longer time as it first computes all the protein ratios to shift spectrum ratios to 1:1. To exemplify this procedure, we only take rat and mouse CERU proteins from `ibspiked_set1`, see figure 3. The resultant noise model is a rough approximation only because of the very limited data, see Breitwieser et al. Supporting Information, submitted, for a real example.

```
R> ib.ceru <- subsetIBSpectra(ibspiked_set1,protein=ceru.proteins,
                             direction="exclude others",
                             specificity="reporter-specific")
R> nm.ceru <- NoiseModel(ib.ceru,one.to.one=FALSE,pool=TRUE)
```

```
3 proteins with more than 10 spectra, taking top 50.
[1] 1.000e-10 4.474e-01 2.057e-01
```

3.4 Protein and peptide ratio calculation

`estimateRatio` calculates the relative abundance of a peptide or protein in one tag compared to another. It calculates a weighted average (after outlier removal) of the spectrum ratios. The weights are the inverse of the spectrum ratio variances. It requires a `IBSpectra` and `NoiseModel` object and definitions of `channel1`, `channel2`, and the protein or peptide. The result is `channel2/channel1`.

```
R> ## Calculate ratio based on all spectra of peptides specific
R> ## to CERU_HUMAN, CERU_RAT or CERU_MOUSE. Returns a named
R> ## numeric vector.
R> 10^estimateRatio(ibspiked_set1,noise.model,
                   channel1="114",channel2="115",
                   protein=ceru.proteins)['lratio']

lratio
0.9272

R> ## If argument 'combine=FALSE', estimateRatio returns a data.frame
R> ## with one row per protein
R> 10^estimateRatio(ibspiked_set1,noise.model,
                   channel1="114",channel2="115",
                   protein=ceru.proteins,combine=FALSE)[,'lratio']

P00450 P13635 Q61147
1.0564 1.8324 0.5067

R> ## spiked material channel 115 vs 114:
R> ## CERU_HUMAN (P00450): 1:1
R> ## CERU_RAT (P13635): 2:1 = 2
R> ## CERU_MOUSE (Q61147): 5:10 = 0.5
```

```
R> png("fig_maplot_noisemodel.png",width=3.3,height=3,units="in",res=600,points=8)
R> maplot(ib.background,noise.model=c(noise.model,nm.ceru),
          channel1="114",channel2="115",ylim=c(0.2,5),
          main="95% CI noise model")
R> dev.off()
```

pdf
2




Figure 3: Red lines denote the 95 % confidence interval as estimated by the noise model on background proteins. The blue line is estimated as non 1:1 noise model based on only spectra of CERU proteins.


```

R>
R> ## Peptides shared between rat and mouse
R> pep.shared <- peptides(proteinGroup(ibspiked_set1),
                          c(ceru.rat,ceru.mouse),set="intersect",
                          columns=c('peptide','n.shared.groups'))
R> ## remove those which are shared with other proteins
R> pep.shared <- pep.shared$peptide[pep.shared$n.shared.groups==2]
R> ## calculate ratio: it is between the rat and mouse ratios
R> 10^estimateRatio(ibspiked_set1,noise.model,
                   channel1="114",channel2="115",
                   peptide=pep.shared)['lratio']

lratio
0.6291

```

When examining the global differences and differences in between classes, `proteinRatios` can be used. It is also suitable to inspect sample variability. The argument `cl` can be used to define class labels. If `method='interclass'` or `intraclass` and `summarize=TRUE`, `proteinRatios` return a single summarized ratio across and within classes, resp..

```

R> protein.ratios <- proteinRatios(ibspiked_set1,noise.model)
R> str(protein.ratios)

'data.frame':      966 obs. of  9 variables:
 $ lratio      : num  -0.0509 -0.0163 -0.0154 0.0344 0.0339 ...
 $ variance    : num  0.000903 0.000671 0.000619 0.000646 0.000534 ...
 $ n.spectra   : num  189 189 189 189 189 189 5 5 5 6 ...
 $ p.value.rat : num  0.0452 0.2646 0.2679 0.0879 0.0712 ...
 $ p.value.sample: num  NA NA NA NA NA NA NA NA NA NA ...
 $ is.significant: num  NA NA NA NA NA NA NA NA NA NA ...
 $ protein     : chr   "136429" "136429" "136429" "136429" ...
 $ r1          : chr   "114" "114" "114" "115" ...
 $ r2          : chr   "115" "116" "117" "116" ...
 - attr(*, "combn.method")= chr "global"
 - attr(*, "symmetry")= logi FALSE
 - attr(*, "sign.level.rat")= num 0.05
 - attr(*, "sign.level.sample")= num 0.05
 - attr(*, "variance.function")= chr "maxi"
 - attr(*, "combine")= logi FALSE
 - attr(*, "reverse")= logi FALSE

R> ## defined class 114 and 115 as class 'T', 116 and 117 as class 'C'
R> classLabels(ibspiked_set1) <- c("T","T","C","C")
R> proteinRatios(ibspiked_set1,noise.model,protein=ceru.proteins,
                cl=classLabels(ibspiked_set1),method="interclass",
                summarize=T)[,c("protein","lratio","variance")]

```

	protein	lratio	variance
1	P00450	0.009653	0.0004753
2	P13635	0.602924	0.0508894
3	Q61147	-0.559911	0.0477306

3.5 Protein ratio distribution and selection

To examine differentially expressed proteins, we usually use both sample variability information (random protein ratios) as a *fold-change* constraint, and ratio variance as a *precision* constraint. For an experimental setup with a different class in each reporter tag, a sample distribution should be learned in advance by generating all possible ratio pairs between samples of the same class. A Cauchy distribution fits accurately this type of random protein ratio distribution: Cauchy is displayed in red, Gaussian in blue. In the case of `ibspiked_set1`, the many 1:1 proteins provide us with adequate data to learn the random protein ratio distribution.

```
R> #protein.ratios <- proteinRatios(ibspiked_set1,noise.model)
R> protein.ratiodistr.wn <- fitWeightedNorm(protein.ratios[, 'lratio'],
                                           weights=1/protein.ratios[, 'variance'])
R> protein.ratiodistr.cauchy <- fitCauchy(protein.ratios[, "lratio"])
```

Now, when supplying a `ratiodistr` parameter to `estimateRatio` and `proteinRatios`, sample and signal p-values are calculated, what we illustrate in the code below

```
R> rat.list <-
  estimateRatio(ibspiked_set1,noise.model=noise.model,channel1="114",channel2="115",
               protein=reporterProteins(proteinGroup(ibspiked_set1)),combine=F,
               ratiodistr=protein.ratiodistr.cauchy)
R> rat.list[rat.list[, "is.significant"]==1,]
```

	lratio	variance	n.spectra	p.value.rat	p.value.sample	is.significant
P13635	0.2630	0.0063314	249	5.261e-04	0.02207	1
Q61147	-0.2953	0.0009033	157	2.049e-23	0.01934	1
<NA>	NA	NA	NA	NA	NA	NA
<NA>	NA	NA	NA	NA	NA	NA
<NA>	NA	NA	NA	NA	NA	NA

3.6 Detection of proteins with no specific peptides

It is well known that MS analysis only reveals the presence of so-called protein groups, defined as sets of proteins identified by the same set of peptides. The protein that contains all the peptides is the group reporter (there are possibly several group reporters) and if it has one specific peptide at least then its presence in the sample is certain. The status of the other proteins in the group is in general impossible to determine. When quantitative information is available, there is a potential to elucidate the structure of part of the protein groups.

In the example below, a subset `IBSpectra` object is created, containing only peptides shared between `CERU_RAT` and `CERU_MOUSE`, and those specific to `CERU_RAT`.

```

R> limits=seq(from=-0.5,to=0.5,by=0.001)
R> curve.wn <- data.frame(x=limits,y=d(protein.ratiodistr.wn)(limits))
R> curve.cauchy<-data.frame(x=limits,y=d(protein.ratiodistr.cauchy)(limits))
R> g <- ggplot(data.frame(protein.ratios),aes(x=lratio)) +
  geom_histogram(colour = "darkgreen", fill = "white",aes(y=..density..),
    binwidth=0.02) + geom_rug() +
  geom_line(data=curve.wn,aes(x=x,y=y),colour="blue") +
  geom_line(data=curve.cauchy,aes(x=x,y=y),colour="red")
R> print(g)

```




Figure 4: Histogram of all protein ratios in `ibspiked_set1`. A fit with a Gaussian and Cauchy probability density function is shown in blue and red, respectively.

```

R> ## peptides shared between CERU_RAT and CERU_MOUSE have been computed before
R> pep.shared

[1] "AGLQAFFQVR"      "DNEEFLESNK"      "DTANLFPHK"      "EMGPTYADPVCLSK"
[5] "ETFTYEWTVPK"     "GSLLADGR"        "KGSLLADGR"      "LYHSHVDAPK"
[9] "NMATRPYSLHAHGVK" "RDTANLFPHK"      "VFFEQGATR"

R> ## peptides specific to CERU_RAT
R> pep.rat <- peptides(proteinGroup(ibspiked_set1),protein=ceru.rat,
                      specificity="reporter-specific")
R> ## create an IBSpectra object with only CERU_RAT and shared peptides
R> ib.subset <- subsetIBSpectra(ibspiked_set1,
                              peptide=c(pep.rat,pep.shared),direction="include")
R> ## calculate shared ratios
R> sr <- shared.ratios(ib.subset,noise.model,
                      channel1="114",channel2="117",
                      ratiodistr=protein.ratiodistr.cauchy)
R> sr

      reporter.protein protein2 ratio1 ratio1.var n.spectra.1   ratio2
lratio      P13635   Q61147 0.9538    0.0117         250 -0.002762
      ratio2.var n.spectra.2
lratio   0.001242         296
R>

```

4 Report generation

Analysis of reports can be analyzed with **isobar** by usage of the Rscript **create_reports.R**. This script reads command line options and a **properties.conf** file to allow a complete analysis generating \LaTeX (and thus PDF) and Excel reports containing all protein ratios and quality control.

The **properties.conf** file in the current folder overwrite the settings set in the global file in the installation directory. See appendix B for the available settings.

```

R> ## execute to find the path and file location in your installation.
R> system.file("report",package="isobar") ## path
R> list.files(system.file("report",package="isobar")) ## files

```

create_reports.R R script which can be used to create QC and PDF reports It initializes the environment, reads properties and calls **Sweave** on QC and DA Sweave files. Additionally it generates a Excel data analysis report by calling **tab2xls.pl**.

isobar-qc.Rnw Sweave file with quality control plots.

```
R> ## plot significantly different protein groups where 90% CI does not overlap
R> ## CERU_MOUSE and CERU_RAT is detected, as expected.
R> shared.ratios.sign(sr,z.shared=1.282)
```

	reporter.protein	protein2	n.spectra.1	n.spectra.2	proteins
1.1	P13635	Q61147	250	296	P13635 \nvs Q61147
1.2	P13635	Q61147	250	296	P13635 \nvs Q61147

	g	ratio	var	n.spectra	id
1.1	reporter	0.953810	0.011703	> 10	1
1.2	member	-0.002762	0.001242	> 10	1




Figure 5: Peptides of spiked ceruplasmins have significantly different ratios between groups. Group *reporter* consists of peptides specific to CERU_RAT (P13635), group *member* are peptides shared between CERU_RAT and CERU_MOUSE (Q61147).

isobar-analysis.Rnw Sweave file for generating a data analysis report with the list of all protein ratios and list of significantly different proteins.

properties.conf Default configuration for `create_reports.R`. It is parsed as R code.

report-utils.R Helper R functions used in Sweave documents.

report-utils.tex Helper L^AT_EX functions used in Sweave documents.

A File formats

A.1 ID CSV file format

The Perl parsers create ID CSV files - identification information for all matched spectra without quantitative information. You can create your own parser, the resulting file should be tab-delimited and contain the following columns. Only bold columns are obligatory. The information is redundant - that means if a peptide may stem from two different proteins the information of the identification is repeated.

accession	Protein AC
peptide	Peptide sequence
modif	Peptide modification string
charge	Charge state
theo.mass	Theoretical peptide mass
exp.mass	Experimentally observed mass
parent.intens	Parent intensity
retention.time	Retention time
spectrum	Spectrum identifier
search.engine	Protein search engine and score

A.2 IBSpectra CSV file format

IBSpectra file format has the same columns as the ID CSV format and additionally columns containing the quantitation information, namely *Xtagname_mass* and *Xtagname_ions*, for mass and intensity of each tag *tagname*. Below an example of the further columns for an iTRAQ 4plex IBSpectra.

X114_mass	reporter ion mass
X115_mass	reporter ion mass
X116_mass	reporter ion mass
X117_mass	reporter ion mass
X114_ions	reporter ion intensity
X115_ions	reporter ion intensity
X116_ions	reporter ion intensity
X117_ions	reporter ion intensity

B properties.conf for report generation

```
R> cat(readLines(system.file("report","properties.conf",package="isobar")),sep="\n")

##
## Isobar properties.conf file
##   for automatic report generation
##
## It is standard R code and parsed using system.file

## Isobaric tagging type. Use one of the following:
# type='iTRAQ4plexSpectra'
# type='iTRAQ8plexSpectra'
# type='TMT2plexSpectra'
# type='TMT6plexSpectra'
type=NULL
isotope.impurities=NULL

## Name of project, by default the name of working directory
## Will be title and author of the analysis reports.
name=basename(getwd())
author="isobar R package"
```

```

ibspectra=paste(name,"ibspectra.csv",sep=".")

## Where should cached files be saved?
# cachedir="cache"
cachedir="."

## An ibspectra object can be generated from peaklists and identifications.

## peaklist files, by default all mgf file in directory
peaklist=list.files(pattern="*\\.mgf")
## id files, by default all id.csv files in directory
identifications=list.files(pattern="*\\.id.csv")
## mapping files, for data quantified and identified with different but
## corresponding spectra. For example corresponding HCD-CID files.

## masses and intensities which are outside of the 'true' tag mass
## +/- fragment.precision/2 are discarded
fragment.precision=0.01
## filter mass outliers
fragment.outlier.probab=0.001

readIBSpectra.args = list(
  mapping.file=NULL
)

correct.isotope.impurities=TRUE

normalize=TRUE
normalize.use.protein=NULL
normalize.exclude.protein=NULL
normalize.function=median
normalize.exclude.set = list (seppro_igy14=c(
  "P02763", # Alpha1-Acid Glycoprotein
  "P01009-1", # Alpha1-Antitrypsin
  "P19652", # Alpha1-Acid Glycoprotein
  "P01023", # Alpha2-Macroglobulin
  "P02768-1", # Albumin
  "P02647", # HDL: Apolipoprotein A1
  "P02652", # HDL: Apolipoprotein A1
  "P04114", # LDL: Apolipoprotein B
  "P01024", # Complement C3
  "P02671-1", # Fibrinogen
  "P00738", # Haptoglobin
  "P01876", # IgA 1
  "P01877", # IgA 2
  "P01857", # IgG 1
  "P01859", # IgG 2
  "P01860", # IgG 3
  "P01861", # IgG 4
  "P01871-1", # IgM
  "P02787" # Transferrin

```

```

));

use.na=FALSE

## the parameter noise.model can be either a NoiseModel object or a file name
data(noise.model.hcd)
noise.model=noise.model.hcd
## If it is a file name, a noise model is estimated as non one-to-one and saved
## into the file. otherwise, the noise model is loaded from the file
# noise.model="noise.model.rda"

## Certain channels can be defined for creation of a noise model
## e.g. if the first and second channel are technical repeats
## If NULL, all channel combinations are taken into account when creating a
## noise model.
noise.model.channels=NULL
noise.model.minspectra=50

summarize=FALSE
combn.method="interclass"
## class labels. Must be of type character and of same length as number of channels
## I. e. 4 for iTRAQ 4plex, 6 for TMT 6plex
## Example for iTRAQ 4plex:
# class.labels=as.character(c(1,0,0,0))
# class.labels=c("Treatment","Treatment","Control","Control")
class.labels=NULL
combn=NULL

## Analysis report sections: Significant proteins and protein details
show.significant.proteins=FALSE
show.protein.details=TRUE

ratios.opts = list(
  sign.level.sample=0.01,
  sign.level.rat=0.01)

quant.w.groupeptides=c("bcrabl","bcrabl,bcrabl_p185,bcrabl_t315i")

min.detect=NULL

database="Uniprot"
preselected=c()

ratiodistr.summarize=FALSE
ratiodistr.summarize.method="global"

write.qc.report=TRUE
write.report=TRUE
write.xls.report=TRUE

```



```
sum.intensities=FALSE
regen=FALSE

scratch=list()
```

C Dependencies

C.1 L^AT_EX and PGF/TikZ

L^AT_EX is a high-quality typesetting system; it includes features designed for the production of technical and scientific documentation. It is available as free software¹. PGF is a T_EX macro package for generating graphics. It comes with a user-friendly syntax layer called TikZ².

L^AT_EX is used for creating PDF analysis reports, with the PGF package creating the graphics. Go to <http://www.latex-project.org> to get information on how to download and install a L^AT_EX system and packages.

C.2 Perl

Perl is a high-level, general-purpose, interpreted, dynamic programming language. Perl is required for two tasks:

- Conversion of Pidres XML and Mascot DAT files to ID CSV format;
- Creation of Microsoft Excel format data analysis report.

Go to <http://www.perl.org> to download and get help on the installation of Perl on your Operating System. For file format conversion, perl module `Statistics::Lite` is required. For Excel export `Spreadsheet::WriteExcel`. All Perl scripts are in the subdirectory `pl` of the *isobar* package installation.

```
R> ## execute to find the path and file location in your installation.
R> system.file("pl",package="isobar") ## path
R> list.files(system.file("pl",package="isobar")) ## files
```

`mascotParser2.pl` and `pidresParser2.pl` convert from respective protein search output-files to a XML file format, which can be converted into a CSV file readable by *isobar* by using `psx2tab2.pl`.

`mascotParser2.pl` converts from Mascot format, and requires the file `modifconv.csv` as a definition of modification names. `pidresParser2.pl` converts from Phenyx output and requires the file `parsersConfig.xml`. `tab2xls.pl` converts csv file to different sheets of an Excel spreadsheet.

```
R> ## execute on your system
R> system(paste("perl",system.file("pl","mascotParser2.pl",package="isobar"),
  "--help"))
R> print(paste("perl",system.file("pl","pidresParser2.pl",package="isobar"),
  "--help"))
```

¹<http://www.latex-project.org>

²<http://sourceforge.net/projects/pgf>

D Session Information

The version number of R and packages loaded for generating the vignette were:

```
R> toLatex(sessionInfo())
```

- R version 2.13.1 (2011-07-08), x86_64-unknown-linux-gnu
- Locale: LC_CTYPE=en_US.UTF-8, LC_NUMERIC=C, LC_TIME=en_US.UTF-8, LC_COLLATE=en_US.UTF-8, LC_MONETARY=C, LC_MESSAGES=en_US.UTF-8, LC_PAPER=en_US.UTF-8, LC_NAME=C, LC_ADDRESS=C, LC_TELEPHONE=C, LC_MEASUREMENT=en_US.UTF-8, LC_IDENTIFICATION=C
- Base packages: base, datasets, graphics, grDevices, grid, methods, stats, utils
- Other packages: Biobase 2.12.2, distr 2.3.2, ggplot2 0.8.9, isobar 0.2.5, plyr 1.5, proto 0.3-9.1, reshape 0.8.4, sfsmisc 1.0-14, startupmsg 0.7.1, SweaveListingUtils 0.5
- Loaded via a namespace (and not attached): biomaRt 2.6.0, digest 0.4.2, RCurl 1.5-0, tools 2.13.1, XML 3.2-0